

AeroDawg

Canine Aerosol Chamber

VETERINARY

LIGHT WEIGHT, METERED DOSE INHALER (MDI) COUPLING CHAMBER FOR DOGS

Chronic bronchitis is one of the most common respiratory diseases in dogs, but is difficult to diagnose as it can mimic other pathologies. Dogs which experience chronic inflammation are often sensitive to environmental and household allergens.

Chronic bronchitis in dogs is caused by inflammation of the airways, resulting in chronic cough and exercise intolerance. Common signs of chronic bronchitis include;

- Daily Coughing**
- Gagging or swallowing after coughing**
- Difficulty breathing**
- Wheezing for a duration of 2-3 months or longer**



Features

- Target specific lung delivery of medication.
- Controls and alleviates airway inflammation and airway obstruction.
- Patented 'Flow-Vu' indicator enables accurate counting of inhalations.
- Exhalation valve directs the flow of air away from your dog to prevent irritation.
- Two mask sizes to accommodate all breeds and sizes of dogs.
- Enables multiple medications to be administered via the one device.

Benefits

- Accurate directing of aerosolised pharmaceuticals to the lungs and small airways.
- Improved drug safety and efficacy by reducing total therapeutic dose required.
- Convenient, portable chamber.
- Simple to use and easy to clean.
- Reduces the risk of side effects and stresses associated with pills, syrups or injections.

DIRECTIONS FOR USE/CONTRAINDICATIONS:

Able to be used as the coupling device for aerolised Metered Dose Inhalers (MDI's). By depressing the MDI, aerolised particles are delivered into the chamber and inhaled. Gently apply the mask to the muzzle and upon inhalation, the 'Flow-Vu' indicator will flutter toward the dog signalling inhalation.

CONSTRUCTION:

Unique, durable, anti-static polymer chamber with flexible nostril cup.

Ranvet Pty Ltd

10-12 Green Street, Banksmeadow NSW, 2019 Australia

Ph: 612 9666 1744 . Fax: 612 9666 1755

Web: www.ranvet.com.au . Email: info@ranvet.com.au