

AeroKat

Feline Aerosol Chamber

VETERINARY

LIGHT WEIGHT, METERED DOSE INHALER (MDI) COUPLING CHAMBER FOR CATS

Approximately **1 in 100** cats have asthma. Additionally, symptoms may often be confused with hairballs. Asthma usually affects cats between 2-8 years of age with Siamese and Himalayan breeds having a higher incidence of respiratory conditions.

If your cat experiences any of the following symptoms they may be suffering from feline asthma;

- Frequent wheezing.**
- Persistent coughing.**
- Open mouthed or laboured breathing.**



Features

- Target specific lung delivery of medication.
- Controls and alleviates airway inflammation and airway obstruction.
- Patented 'Flow-Vu' indicator enables accurate counting of inhalations.
- Exhalation valve directs the flow of air away from your cat to prevent irritation.
- Two mask sizes to accommodate all breeds and sizes of cats.
- Enables multiple medications to be administered via the one device.

Benefits

- Accurate directing of aerosolised pharmaceuticals to the lungs and small airways.
- Improved drug safety and efficacy by reducing total therapeutic dose required.
- Convenient, portable chamber.
- Simple to use and easy to clean.
- Reduces the risk of side effects and stresses associated with pills, syrups or injections.

DIRECTIONS FOR USE/CONTRAINDICATIONS:

Able to be used as the coupling device for aerosolised Metered Dose Inhalers (MDI's). By depressing the MDI, aerosolised particles are delivered into the chamber and inhaled. Gently apply the mask to the muzzle and upon inhalation, the 'Flow-Vu' indicator will flutter toward the cat signalling inhalation.

CONSTRUCTION:

Unique, durable, anti-static polymer chamber with flexible nostril cup.

Ranvet Pty Ltd

10-12 Green Street, Banksmeadow NSW, 2019 Australia

Ph: 612 9666 1744 . Fax: 612 9666 1755

Web: www.ranvet.com.au . Email: info@ranvet.com.au