

AeroDawgTM*



VETERINARY

LIGHT WEIGHT, METRED DOSE INHALER (MDI) COUPLING CHAMBER FOR DOGS.

Chronic bronchitis is one of the most common respiratory diseases in dogs, but is difficult to diagnose as it can mimic other pathologies. Dogs which experience chronic inflammation are often sensitive to environmental and household allergens.

Chronic bronchitis in dogs is caused by inflammation of the airways, resulting in chronic cough and exercise intolerance. Common signs of chronic bronchitis include;

- Daily coughing.**
- Gagging or swallowing after coughing.**
- Difficulty breathing.**
- Wheezing for a duration of 2-3 months or longer.**

Features

- Target specific lung delivery of medication.
- Controls and alleviates airway inflammation and airway obstruction.
- Patented 'Flow-Vu' indicator enables accurate counting of inhalations.
- Exhalation valve directs the flow of air away from your dog to prevent irritation.
- Two mask sizes to accommodate all breeds and sizes of dogs.
- Enables multiple medications to be administered via the one device.

Benefits

- Accurate directing of aerosolised pharmaceuticals to the lungs and small airways.
- Improved drug safety and efficacy by reducing total therapeutic dose required.
- Convenient, portable chamber.
- Simple to use and easy to clean.
- Reduces the risk of side effects and stresses associated with pills, syrups or injections.

DIRECTIONS FOR USE/CONTRAINDICATIONS

Able to be used as the coupling device for aerosolised Metered Dose Inhalers (MDI's). By depressing the MDI, aerosolised particles are delivered into the chamber and inhaled. Gently apply the mask to the muzzle and upon inhalation, the 'Flow-Vu' indicator will flutter toward the dog signalling inhalation.

CONSTRUCTION

Unique, durable, anti-static polymer chamber with flexible nostril cup.